






Objectives:

This experiment aims to analyze the composition, quality, and relative compositions of canola oil (CO), palm oil (PO), and a 2:1 CO-PO mixture in relation to the molecular structures of the fatty acids present in both oils via ^1H -NMR and FTIR.

Reagents:

<i>Reagent Name and Structure</i>	<i>Physical and Chemical Properties</i>	<i>Safety Precautions</i>
Deionized Water (H_2O) 	<ul style="list-style-type: none"> - colorless, odorless liquid - MW = 18.02 g/mol 	Non-irritant, generally safe to handle Disposal: <u>sink</u>
Canola Oil (CO)	<ul style="list-style-type: none"> - light yellow liquid with a slight oil odor - less dense than water (specific gravity ~ 0.9) - comprised of oleic, linolenic, and linoleic acids (91%) 	<ul style="list-style-type: none"> - somewhat flammable - wash with water and soap upon direct contact - keep away from drains Disposal: <u>I104</u>
Palm Oil (PO)	<ul style="list-style-type: none"> - light yellow solid with a characteristic, slight odor - has melting point range from 25 to 40 $^{\circ}\text{C}$ - less dense than water (specific gravity ~ 0.9) - comprised of palmitic acid (44%) and oleic and linolenic acids (50%) 	<ul style="list-style-type: none"> - somewhat flammable - wash with water and soap upon direct contact - keep away from drains Disposal: <u>I104</u>
Chloroform-d (CDCl_3) 	<ul style="list-style-type: none"> - clear, odorless liquid - somewhat soluble in water - denser than water (density = 1.500 g/cm3) - does not give ^1H-NMR signals when used as a solvent 	<ul style="list-style-type: none"> - can cause serious eye and skin irritation - toxic if inhaled - can cause damage to organs (liver and kidney) if ingested - can cause cancer / carcinogenic - Disposal: <u>G703</u>

<p>Tetramethylsilane (TMS)</p> 	<ul style="list-style-type: none"> - colorless, odorless liquid - low boiling point from 26 to 28 °C - less dense than water (density = 0.648 g/cm³) - yields reference signal at $\delta = 0$ ppm when mixed with sample for ¹H-NMR 	<ul style="list-style-type: none"> - extremely flammable in both liquid and vapor forms - Disposal Code: <u>G704</u>
<p>Potassium Bromide (KBr)</p> 	<ul style="list-style-type: none"> - appears as a white, odorless, crystalline solid - very soluble in water; KBr crystals are destroyed in the presence of moisture - transparent to IR radiation; used as the material for sample holder in IR spectroscopy 	<ul style="list-style-type: none"> - causes skin and eye irritation upon direct contact - should be handled with gloves on - can cause respiratory irritation if KBr dust is inhaled - harms aquatic life - Disposal Code: <u>D499</u>
<p>Acetone</p> 	<ul style="list-style-type: none"> - colorless liquid with a mildly sweet, fruity odor - soluble in water but will not dissolve KBr crystals - may be used to clean IR-related equipment (ask technicians for cleaning steps) 	<ul style="list-style-type: none"> - highly flammable in liquid and vapor form - causes serious eye irritation - causes respiratory irritation if inhaled - Disposal Code: <u>G704</u>

Required Equipment:

- IRAffinity Machine
- Bruker NMR Machine
- Glassware for Producing 2:1 CO-PO Mixture (e.g. graduated cylinder)
- Glassware for Loading Oil Samples, TMS, and CDCl₃ for NMR Spectroscopy
- Sandwich and Pellet Setups for Liquid and Solid Samples for IR Spectroscopy

Literature Data

Table 1. FTIR Absorptions Present in Edible Fats and Oils from $\nu = 4000$ to 600 cm^{-1} (Ye & Meng, 2022)

#	Bond(s) Involved	Wavenumber ν (cm^{-1})	Type of Vibration
1	C-H	720	C-H stretching (rocking)
2	C-O	1166	C-O stretching
3	C-H	1454	C-H bending
4	C=O	1749	C=O stretching
5	C-H	2856	C-H symmetric stretching
6	Alkene C-H	2937	Alkene C-H stretching

Table 2. ^1H -NMR Chemical Shifts of Signals Elicited by Edible Oils in CDCl_3 (Siudem et al., 2022)

#	Compounds	Group Involved	Chemical Shift δ (ppm)
1	All fatty acids except linolenic acid	$-\text{CH}_2-\text{CH}_3$	0.89
2	Linolenic acid	$-\text{CH}=\text{CH}-\text{CH}_2-\text{CH}_3$	0.97
3	Fatty acids with acyl chains	$-(\text{CH}_2)_n-$	1.31
4	Fatty acids with acyl chains	$-\text{CH}_2-\text{CH}_2-\text{COOH}$	1.61
5	Mono- and polyunsaturated fatty acids	$-\text{CH}_2-\text{CH}=\text{CH}-$	2.04
6	Unsaturated fatty acids with acyl chains	$-\text{CH}_2-\text{COOH}$	2.31
7	Linoleic and linolenic acids	$-\text{CH}=\text{CH}-\text{CH}_2-\text{CH}=\text{CH}-$	2.77
8	Triacylglycerols	$-\text{CH}_2-\text{OCO}-$	4.27
9	Triacylglycerols	$-\text{CH}-\text{OCO}-$	5.28
10	Mono- and polyunsaturated fatty acids	$-\text{CH}=\text{CH}-$	5.35

Table 3. Relative Compositions of Canola and Palm Oils (Fatty Acid Composition, 2018)

Fatty Acid	Lipid Number	wt% in Palm Oil	wt% in Canola Oil (Rapeseed Oil)
Lauric	12:0	0.3	-
Myristic	14:0	1.0	-
Palmitic	16:0	44	4
Stearic	18:0	4	1.5
Oleic	18:1	40	60
Linoleic	18:2	10	20
Linolenic	18:3	0.2	11
Arachidic	20:0	-	0.5
Gadoleic Eikosadienic	20:1+2	-	1.5
Behenic	22:0	-	0.4
Erucic	22:1	-	0.4

General Observations:

References:

- Fatty Acid Composition of Selected Oils and Fats.* (2018, June 22). Fediol. <https://www.fediol.eu/data/fatty%20acids.pdf>.
- Siudem, P., Zielinska, A., Paradowska, K. (2022). Application of ¹H NMR in the study of fatty acids composition of vegetable oils. *Journal of Pharmaceutical and Biomedical Analysis*, 212, 114568. <https://doi.org/10.1016/j.jpba.2022.114658>
- Ye, Q., Meng, X. (2022). Highly efficient authentication of edible oils by FTIR spectroscopy coupled with chemometrics. *Food Chemistry*, 385, 132661. <https://doi.org/10.1016/j.foodchem.2022.132661>
- chemistryconnection.com/sds/data/pdf/Canola_Oil_SDS_CC.pdf
- chemos.de/import/data/msds/GB_en/8002-75-3-A0050032-GB-en.pdf
- sigmaaldrich.com/PH/en/sds/aldrich/494275?userType=undefined
- sigmaaldrich.com/PH/en/sds/aldrich/87921?userType=undefined
- fishersci.com/content/dam/fishersci/en_US/documents/programs/education/regulatory-documents/sds/chemicals/chemicals-p/S25478.pdf